



Aixplorer® Tips

Evaluation of Liver Fibrosis using ShearWave[™] Elastography on Aixplorer[®]

SSID02892-01 November 2012





From a clinical standpoint, staging liver fibrosis is of major importance^{1, 2}:

- to make a prognosis
- to follow up the evolution of chronic liver diseases (cirrhosis or hepatitis)
- to monitor antifibrotic treatments

The stage of the fibrosis is correlated to the stiffness of the tissue, which is also correlated to the velocity of the shear wave

ShearWaveTM Elastography (SWETM) performed with Aixplorer[®] is a 2D imaging mode^{1,3} that provides a view of the organ and a map of the shear wave velocity over the region of interest. With the Adjustable Numerical Scale, the map can be fine-tuned to highlight the stiffer or the softer parts of the ROI.

Apart from fibrosis, many studies demonstrated the influence of several clinical factors on liver stiffness:

- 1. Respiration⁴, deep breath
- 2. Central venous pressure⁵
- 3. Intrahepatic cholestasis⁶⁻⁷
- 4. Hepatic necro-inflammatory activity (steatosis⁸ for example)
- 5. Peliosis hepatitis (affection of the liver parenchyma vasculature)
- 6. Thrombosis of hepatic vein (clot)
- 7. Congestive hepatopathy⁹...

These should be considered when assessing liver stiffness.

The known limitations of conventional ultrasound examination also apply to the SWE[™] mode: narrow intercostal spaces, thick layer of fat...

Noteworthy: the presence of ascites is not a limitation for the evaluation of liver fibrosis with SWE^{TM} .

Reference

- 1. Quantitative viscoelasticity mapping of human liver using supersonic shear imaging: preliminary in vivo feasibility study. *Muller M, Gennisson JL, Deffieux T, Tanter M, Fink M.* Ultrasound Med Biol. 2009 Feb;35(2):219-29.
- 2. Management of Hepatitis C: June 2002. *National Institute of Health Consensus Development Conference Statement.* Hepatology 2002; 36[5 suppl 1]:S3–S20
- 3. Supersonic shear imaging: A new technique for soft tissues elasticity mapping. *Bercoff J, Tanter M, Fink M.* IEEE Trans Ultrason Ferroelectr Freq Control 2004c;51(4):396–409.
- 4. The effect of the respiratory cycle on liver stiffness values as measured by transient elastography. Yun MH, Seo YS, Kang HS, Lee KG, Kim JH, An H, Yim HJ, Keum B, Jeen YT, Lee HS, Chun HJ, Um SH, Kim CD, Ryu HS. J Viral Hepat. 2010 Oct 4.
- Liver stiffness is directly influenced by central venous pressure. Millonig G, Friedrich S, Adolf S, Fonouni H, Golriz M, Mehrabi A, Stiefel P, Pöschl G, Büchler MW, Seitz HK, Mueller S. J Hepatol. 2010 Feb;52(2):206-10.
- 6. Transient elastography. *Della Valle S, Branchi F, Rigamonti C, Fraquelli M.* Recenti Prog Med. 2009 Apr;100(4):195-201.
- 7. Transient elastography (FibroScan): a new tool in hepatology. *Vergniol J, de Lédinghen V.* Presse Med. 2009 Oct;38(10):1516-25.
- 8. Ultrasound-based Hepatic Elastography Origins, Limitations, and Applications, Eric B. Cohen, MD* and Nezam H. Afdhal, MDw (J Clin Gastroenterol 2010;44:637–645).
- 9. Hepatic congestion plays a role in liver stiffness. *Frulio N, Laumonier H, Balabaud C, Trillaud H, Bioulac-Sage P.* Hepatology. 2009 Nov;50(5):1674-5.

Imaging & Quantification Protocol

1. The patient

- He/she lies in the supine position, with the right arm in maximum abduction to make the right hypochondrium accessible.
- Fasting is mandatory.

Aixplorer® Tips

• Normal breathing is recommended.



2. The probe

- Choose the SC6-1 curved probe, with the preset "Liver" in the "Abdominal" application.
- Scan intercostally (probe parallel to the intercostal space within the space) with sufficient gel and minimize rib shadowing. If SWE signal is weak, unstable, and/or not visible, apply a pressure on the probe to open the rib space (improve the acoustic window) and decrease tissue thickness between the probe and the ribs. Contrary to what has been recommended as a rule for most of the organs, a pressure must be applied to the probe when scanning the liver. The ribs will absorb the pressure and the elasticity of the liver will not be impacted.
- The patient can be scanned subcostally, although it is not recommended. If scanned subscostally, apply the slightest pressure required for visualization in B-Mode.
- Slow or even no movement of the probe is recommended to avoid motion artifact and to stabilize the map.

The SWE[™] acquisition of the right lobe is recommended in priority. However, if it has to be performed on the left lobe, apply the slightest pressure required for visualization in B-Mode. Be aware that the SWE[™] acquisition on the left lobe is more subtle and requires no pressure.



3. SWE[™] and the B-Mode image

The brightness of the B-Mode image is related to the best acoustic window for the SWE velocity map. **Ensure the optimal B-Mode image before engaging the SWE™ Mode and placing the SWE™ box:**

- Enlarge the intercostal space and decrease the thickness of the subcutaneous fatty layer, by using correct patient positioning and pressure on the probe.
- Ensure the optimum contact between the probe and the skin.
- Place the probe parallel to the intercostal window to avoid shadowing from the ribs.

4. The SWE[™] Velocity map

Maintain pressure on the probe.

The SWE[™] default settings have been optimized for the assessment of liver fibrosis. Run the first exam with the default settings. Adjust them only if necessary.

- a. The SWE[™] box
 - Move it onto a vessel-free parenchyma.
 - Place it within a zone of uniform parenchyma as defined by the B-Mode image.
 - Avoid placing the SWE[™] box close to the liver capsule.
 - The most robust acquisition is performed between 3 and 7 cm in depth.





b. SWE[™]Optimization (Res/Std/Pen)

If a lack of SWE[™] signal is observed, turn to "Penetration"

c. Freezing the image

- The patient should hold his/her breath for at least 4 seconds during the expiration phase.
- This delay allows sufficient filling of the SWE box and stabilization of the image in a nomovement context.
- Freeze

Aixplorer® Tips

≥3.2 m/s

≥3.2 m/s

2.8

≥7.1 m/s

≥2.2 m/s

2 (

d. SWE™ Adjustable Numerical Scale (Live or frozen)

Default value: 0-3.2 m/s

- Increasing the scale enables the velocities higher than 3.2 m/s to be spread over a wider range, hence to better discriminate the areas with different velocities higher than the default value.
- Decreasing the scale enables the velocities lower than 3.2 m/s to be spread over a wider range, hence to better discriminate the areas with different velocities lower than the default value.





c. Assessing the reliability of the SWE $^{\scriptscriptstyle\rm TM}$ acquisition

Default ANS: 3.2 m/s

- Repeat this procedure 3 times to acquire 3 valid, independent SWE[™] images of the same scanning view.
- If the 3 acquisitions are not similar, it is recommended to consider the test a failure.





Examples of acquisitions that should be considered unsuccessful

B-mode image

Shadowing, poor image quality, rushed acquisition.

Non valid for the SWETM Velocity map



SWE[™] Velocity map

SWE box underneath the capsule, liver motion due to the patient's respiration, lack of SWE signal, rushed acquisition.

Non valid SWE[™] acquisition



SuperSonic Imagine France

Les Jardins de la Duranne Bat E & F 510, rue René Descartes 13857 Aix-en-Provence Cedex France ① +33 (0)4 42 99 24 32 봄 +33 (0)4 42 52 59 21 ⊠ contactsFR@supersonicimagine.fr Hotline (du lundi au vendredi 8h30-17h30) : 04 26 70 02 70

SuperSonic Imagine, Inc USA

SuperSonic Imagine Asian Distribution Network

Les Jardins de la Duranne Bat F 510, rue René Descartes 13857 Aix-en-Provence, France ① +33 (0)4 42 99 24 32 圖 +33 (0)4 42 52 59 21 🖂 contactsASIA@supersonicimagine.com

www.supersonicimagine.com

SuperSonic Imagine Ltd. UK

18, Upper Walk Virginia Water Surrey GU25 4SN UK ◎ +44 (0)845 643 4516 ⊠ contactsUK@supersonicimagine.com

SuperSonic Imagine GmbH Deutschland

Dietlindenstr. 15 80802 München Deutschland ⊕ +49 89 36036 844 ♣ +49 89 36036 700 ⊠ contactsDE@supersonicimagine.com

